



**Built 4 Me**

**04/2009**





## INTRODUCTION

### **Built 4 Me - as individual as a finger print**

People with a handicap often have completely individual requirements for a wheelchair. We really explore these requirements in great detail. Finding made-to-measure and individual solutions and converting these into reality on our wheelchairs - that is the task of our Built 4 Me department.

Have a look at the Groove range from Quickie, the first powerchairs with modular construction, which fulfil the requirements for design, performance and ease of use, without compromise.

In this brochure you will see some of the special build possibilities in the Groove range. You can either order these parts at the same time that you order your chair, or order them to retrofit afterwards.

Please note that only a few selected special builds are illustrated here.

It is always worth asking our team if you have a particular request.





## CONTENTS

### GROOVE M “SPORT” - E-HOCKEY

Groove E-Hockey in various versions	4
-------------------------------------	---

### VARIOUS SPECIAL BUILDS ON THE SEATING SYSTEM - WHAT GOES WITH WHAT?

Perfect Fit seating system with Comfort backrest	6
Abduction wedge on the Comfort seating system	6
Standard sideguard on the Perfect Fit seating system	7
Comfort seat with moving armrest	7

### INDIVIDUAL KÖNIG SEAT - ELECTRICAL OPTIONS

Special König seats	8
53° seat angle module	9
Backrest adjustment with length compensation	9
Various footboard options	10
Plug-free footrest (electric flip-up footrest)	11
Foot control (HMC)	11
Calf pad for 80° footrest hanger	11

### GENERAL SPECIAL BUILDS

16” rear wheel	12
Immobiliser	12
Switch solutions	12
USB connection	13
12V connection	13
Brackets for respiratory equipment	13
JAY3	13
Over-head chin control	14
Side chin control	14
Height adjustable, swing-away console	15
Protos low movement joystick	15



## GROOVE M "SPORT"

In next to no time, you can change the standard Groove M into a Groove M "Sport". The kit for the rear ram protection including ventilator, the quick-release solution to fix the anti-pitch mechanism and pre-tensioning the suspension can all be carried out in around 30 minutes. For E-Hockey, for example, we recommend a 10 kph or 13 kph Groove M.

In addition, you can specify particular colour requests for seat frames, back bar, footrest and rim. Everything else is standard. Let's Groove!

## E - HOCKEY

### "White" E-Hockey Groove

- 10 kph or 13 kph
- Rear ram protection including ventilator
- Blocked suspension
- Fixing for anti-pitch mechanism (quick-release axle)
- 75° offset footrest, continuous footboard
- White seat frame, backrest bar, footrest
- White Easy Max sideguard



### "White" E-Hockey Groove

- As above with Jay 3 backrest
- White-coated rim



### "Bright lime green" E-Hockey Groove

- 10 kph or 13 kph
- Rear ram protection including ventilator
- Blocked suspension
- Fixing for anti-pitch mechanism (quick-release axle)
- 75° offset footrest, continuous footboard
- Seat frame, back bar, footrest in bright lime green
- Easy Max sideguard in bright lime green
- JAY3 backrest
- Rim coated in bright lime green





#### rear ram protection including ventilator

- Colour-coated - as required



#### "Special" ram bar

- Special ram bar for the Groove M "Sport"



#### "Green" E-Hockey Groove

- Bright lime green
- Footrest, seat frame, sideguard, rim and ram bar, all colour-coated
- Jay 3 backrest
- Ram protection, colour-coated, including ventilator



## SEATING SYSTEMS

Whether your requirements are for a standard seating system, a Perfect Fit seating system, a Comfort seating system or Recaro / König seats, there are various special builds we can offer on each of the different seating systems. Seating posture, support and comfort, to stabilise, support and position the client - the most important thing is for the user to sit well and without problems.

It is always worth sending us your individual requests.

### Comfort backrest on the Perfect Fit seating system

- Comfort backrest with armrest from Perfect Fit seating system (to assist side transfer to and from the wheelchair)
- Abduction wedge on the Perfect Fit seating system
- Calf pad on 80° footrest



### Abduction wedge on the Perfect Fit seating system

- Including mobile telephone case





### Standard sideguard on the Perfect Fit seating system

- Standard sideguard on the Perfect Fit seating system



### Comfort seat with moving armrest

- Moving armrest on the Perfect Fit seating system with the Comfort seat
- Backrest adjustment with length compensation
- Plug-free electric footrest



## SEATING SYSTEMS

### Seating means individuality

#### KÖNIG SEAT AND BACKREST WITH LENGTH COMPENSATION

##### König seat (special dimensions)



- König seat (special dimension, to suit the needs of the customer)
- Brace-style belt
- Cut-out for Jay cushion (also for third party products)



##### "Flat" König seat

- Special König seat, which only adds 7 cm to the seat height





### Electrical 53° seat tilt module

- 53° seat tilt module for optimal pressure relief



### Backrest adjustment with length compensation (Shear Reduction)

- Moves backwards on a rail system.

Can be adjusted individually to suit different body dimensions  
(adjustment using side guide)





## FOOTREST

### Footboard “as for Perfect Fit”

- Footboard for centrally mounted footrest width and depth adjustable



### “Special” footboard

- Footboard manufactured individually to meet customer requirements





### Plug-free footrest

- Flip-up, electrical footrest, can be removed without having to first release a connector



### Foot control

- The wheelchair is driven by the foot and electrical options are controlled in this way too (buddy buttons can also be controlled by the foot)



### Calf pad for 80° footrest

- Calf pad (normally with a footrest which is flip-up or centrally mounted) on 80° footrest hanger





## GENERAL

Individual solutions, to meet customer requirements - as unique as a finger print. Here you will get an overview of what is possible.

Here, the same applies - it is always worth asking us

## GENERAL SPECIAL BUILDS

### 16" Drive wheel on Groove F

- 16" drive wheel on Groove Front Wheel Drive for an increased range, and outdoor use.



### Immobiliser

- A jack plug means that you can immobilise the wheelchair




### Switch solution 1

- Smaller switch solution than the usual Buddy Buttons



### Switch solution 2

- Coloured buttons to control the electrical seat options





### JAY3 backrest

The JAY3 positioning backrest on the Groove. Total Comfort – the user is supported and stabilised by the backrest



### USB connection

*USB Docking Station*



### 12 Volt connection



### Brackets for respiratory equipment

The most varied brackets for the most varied respiratory equipment, built to specification.



## CONTROL OPTIONS

Individual solutions, to meet customer requirements - as unique as a finger print. Here you will get an overview of what is possible.

As always, the same applies - it is always worth asking us whatever your requirement.

## CONTROL

### Chin control electrical, swing-away (over head) on the Comfort seat

Chin control electrical, swing-away over the head, with soft ball. Control of menu, electrical options and emergency stop using Buddy Buttons and strip switch



### Chin control, electrical, swing-away (to the side) on the König seat

Chin control, electrical, swing-away to the side, with soft ball. Control of menu, electrical options and emergency stop using Buddy Buttons





Individually height-adjustable console which can be swung away by 90°

The console can be individually height-adjusted and swung outwards by 90° (operated by push button)



Protos low movement joystick

Extremely sensitive joystick control - minimal effort required



Sunrise Medical GmbH & Co.KG  
Kahlbachring 2-4  
69254 Malsch / Heidelberg  
T +49 (0) 7253 980 0  
F +49 (0) 7253 980 111  
[www.sunrisemedical.de](http://www.sunrisemedical.de)

**INNOVATE.<sup>™</sup>**  
**BUILD FOR THE FUTURE**